

Flex_{MD}



THE MOST
POWERFUL PORTABLE
LED DEVICE IN THE WORLD



DISCOVER WHAT IS POSSIBLE

The Flex MD is the gold standard of professional, portable LED Phototherapy, powered by unique precision Tri-Wave MD technology and manufactured in the UK by Aesthetic Technology Ltd.

Trusted by skin care professionals worldwide, the medically certified Flex MD provides unrivalled, clinically proven results for a wide range of indications and combination applications. The versatile flexible design ensures an effective and hygienic treatment delivery for the face and body.



REJUVENATION ACNE PIGMENTATION REDNESS
PSORIASIS WOUND HEALING PAIN PRE/POST TREATMENT

THE POWER OF LIGHT

MOST CLINICALLY PROVEN WAVELENGTHS



Dermalux Blue 415nm


BLUE 415nm is the most proven wavelength to destroy the bacteria which causes acne and improve skin clarity without irritation. Powerful antibacterial action also has a calming effect to reduce redness and inflammation.

Dermalux Red 633nm

RED 633nm supercharges cell energy to accelerate repair and renewal, boost collagen synthesis, increase hydration and enhance circulation to reverse the visible signs of ageing. Ideal for restoring vitality and glow.

Dermalux NIR 830nm

NEAR INFRARED 830nm targets deeper cells to trigger the wound healing process, calm redness and irritation, reduce pigmentation and promote healthy skin. Combine with Red 633nm for advanced rejuvenation benefits.



MOST POWERFUL

Unique Dermalux ceramic LED technology delivers maximum power and minimal heat to ensure optimum performance and consistent treatment outcomes. Results from independent testing demonstrate the Flex MD to be at least three times more powerful than other portable LED devices available.

MOST PROVEN

The Flex MD is a certified medical device, tested against stringent quality, safety and efficacy standards and independently verified to deliver against medical claims. Featuring the most clinically proven wavelengths backed by science, with Dermalux you will simply deliver better results.


REAL RESULTS PAIN FREE



ACNE



REDNESS



PSORIASIS

WHAT'S INCLUDED



ONLINE
CERTIFIED
TRAINING



ONGOING
EDUCATION



LIFETIME
MARKETING
SUPPORT



24/7 CUSTOMER
SUPPORT



TECHNICAL
SUPPORT



2 YEAR
WARRANTY

TECHNICAL SPECIFICATION

WAVELENGTHS (+/- 2nm):

415NM, 633NM, 830NM

CLASSIFICATION:

CLASS IIA MEDICAL DEVICE

TREATMENT TIME (MINS):

1 TO 30

DOSE RANGE (JOULES):

10 TO 40

INTENSITY (M/W):

5 TO 27

DEVICE DIMENSIONS IN BASE UNIT (MM):

W 520 / D 320 / H 300

LED ARRAY DIMENSIONS (MM):

W 280 / D 790 / H 30

DEVICE WEIGHT (KG):

1.5

Don't just take our word for it...

"We have been able to enhance every treatment that we offer with the Flex MD and the visible effects are immediate"

Dr Shirin Lakani, Elite Aesthetics

"The most powerful portable LED device available and the results from the get go, are astonishing"

As recommended by The Telegraph

"The Flex MD has been a game changer for my business"

*Daphne Kapondoro,
Broadland Aesthetics*

10 REASONS TO CHOOSE DERMALUX

1 UNIQUE CUTTING EDGE TECHNOLOGY

2 MOST CLINICALLY PROVEN TREATMENT

3 MEDICALLY CERTIFIED DEVICES

4 MULTI-AWARD WINNING

5 REAL RESULTS, PAIN FREE

6 WORLD LEADING, TRUSTED BRAND

7 MANUFACTURED IN THE UK

8 CONTINUOUS TRAINING AND EDUCATION

9 LIFETIME MARKETING SUPPORT

10 UNRIVALLED PROFITABILITY



MANUFACTURER OF THE YEAR



Aesthetic Technology Ltd.
Tel: 0845 689 1789 E: info@dermaluxled.com
W: dermaluxled.com

@dermaluxled

