

The Lipo Catalogue









Index

	Page
Fluid Infusion	5
Powered Suction	10
Toomey Harvest	15
Luer Lock Harvest	18
Fat Transfer	21
Compression Garments	27
Ancillary Equipment	33

Contact Us

Triangle Surgical Ltd 6-10 Church Farm Courtyard High Street Chalfont St Giles Buckinghamshire HP8 4QH UK Tel: +44 (0)1494 876990

Fax: +44 (0)1494 876991

email: info@trianglesurgical.uk



Welcome

Thank you for choosing Triangle Surgical, we believe that it is our job to make life easy for you. We do this by listening to your requirements, advising you when you ask for suggestions, ensuring you get a sensible price and making sure that what you order is sent when you want it – usually the next day.

We are a long established family run company with staff who have all been with the company for well over ten years – some with over 25 years involvement. Hence we have a significant combined experience in the sector of Cosmetics, Aesthetics and Reconstruction. We take great pleasure and pride in our team's outstanding reputation for quality of service and in-depth product knowledge.

Our focus is on designing, sourcing and marketing superior products that offer real advantage and advancement to the doctor and nursing team.

We make our own cannula enabling us to provide the widest possible range available and also offer a bespoke, made to measure service. Our CETTA (single use) and PRIMO (reusable) cannula are recognised as the market leader in many countries for both choice and quality.

Within this catalogue we showcase a comprehensive range of products which have a proven track record with medical professionals and patients. We offer quality products for all aspects of liposuction and fat transfer along with our optimum comfort compression garments.

You can be confident in the knowledge that all the products we offer are fully compliant and CE certified. We are on the NHS & MIA Register, fully insured and committed to patient safety.

We are Peppol compliant with e-procurement catalogues available via numerous providers, including Nexus and Virtualstock-the Edge.

Our website: www.trianglesurgical.uk













Disposable - Black

The outstanding quality of the CETTA disposable product range is immediately apparent; high grade medical stainless steel, reassuring weight, excellent balance and superb finish. We strongly believe CETTA is not only the best disposable product available worldwide, but is also far superior to the majority of re-usable cannula available.



Cetta handle: POM

CETTA Cannula - identified by black handles.

Certain styles of handles are made using POM - Polyoxymethylene (An exceptional and recyclable engineering thermoplastic).



Cetta Toomey - hub Aluminium, shaft stainless steel



Cetta handle POM, shaft stainless steel



Cetta luer lock hub and shaft stainless steel

CETTA is single packed for your convenience - and cash flow

We will create custom packs for individual surgeon's requirements.

Unlike other disposable brands, the CETTA luer lock hub is metal meaning the quality is exactly the same as the best available re-usable cannula



Ask about our custom packs today





Autoclavable - Gold

Our re-usable Primo range is as versatile as you need it to be, we will create any bespoke style and size upon request at standard list prices. The range shown in this catalogue is a guide to the most popular sizes and styles we currently provide. Many are regular stock items facilitating fast delivery. Made of durable high quality medical grade stainless steel with reassuringly weighty and balanced Aluminium handles available.



Primo handle: Aluminium







Primo handle Aluminium, shaft stainless steel



Feeling Creative?

We often work with surgeons to create their ideal size and tip design of cannula for a given purpose.



We can create bespoke CETTA and Primo from a basic drawing and description

Ask about our bespoke design service today



Fluid Infusion

quality. We are extremely proud to be working with them.

We went straight to the inventor of the Technique!

When It comes to the Tumescent technique we simply could do no better than to partner Dr Jeffrey Klein and his company HK surgical USA.

Often referred to as the Klein technique the use of Tumescent fluid to improve the comfort and recovery of liposuction surgery has become the worldwide standard since Dr Klein's original work in 1986.

In order to provide the very best equipment for his requirements Dr Klein founded HK surgical, his company has become a recognised world leader in the specialist supply of the devices and instrumentation needed for the procedure. We have worked with many devices over the years for both infusion and suction and can confirm that the HK products are exceptional in design and build





Dr Jeffrey Klein

The Complete Tumescent Solution

HK Surgical is the recognised leader in the development of highly specialized Tumescent Local Anesthesia products in the medical device industry.

CAN'T FIND WHAT YOU'RE LOOKING FOR?

WE'RE HAPPY TO HELP!

Call us today on 01494 876990 or email info@trianglesurgical.uk



Infiltration Pumps

Infiltration of large volumes of tumescent anesthesia can optimise patient comfort and procedure result.

However, this can be a laborious job using a hand held syringe and on large cases a pump for infiltration is an excellent solution.

A practical way of greatly reducing expensive theatre time, and reducing the very real physical stresses experienced by the surgeon during a large procedure.

The Klein KIP Pump

The Klein KIP11 pump, although small and manageable, is a little power-house with a 0-600rpm motor producing flow rates of up to 1000ml per minute.

It's a well-considered piece of engineering design with clever touches such as two footswitch ports allowing that one pedal can be positioned either side of the O.R. table for ease and time saving.

The only infiltration pump on the market with a nonconductive roller head is an important additional safety feature.

We think this ultra reliable little device is a very wise investment.

ole, is a wrates

Please get in contact for full spec details.

Product Code

THK-KIP11



The Klein Touch Pump

The new Klein Touch Pump features an LCD touch screen, precision .01 rate increments and the flow rate of up to 1000ml/min. This new well engineered Klein touch pump can be operated from the LCD screen or from one of the three foot pedal ports. It can be pre-programmed and has robust switches with solid and reliable foot pedals.

We strongly believe that this is an efficient, essential, quality device designed to provide consistent results.

Maintenance free brushless motor.

Can be battery operated (optional and can be purchased separately).

Please get in contact for full spec details.

Product Code

THK-KTP



Infusion Needles

High grade stainless steel. 6 spaced holes in a spiral around the tip.

For fluid infusion using luer lock syringes, handles or giving set tubing.



CETTA RANGE	PRIMO RANGE		
Product Code	Product Code	Diameter (mm)	Length (cm)
CETTA INF/N 215	PRIMO INF/N 215	2	15
CETTA INF/N 223	PRIMO INF/N 223	2	23
CETTA INF/N 230	PRIMO INF/N 230	2	30
CETTA INF/N 315	PRIMO INF/N 315	3	15
CETTA INF/N 323	PRIMO INF/N 323	3	23
CETTA INF/N 330	PRIMO INF/N 330	3	30



Product Code	Product Code	Diameter (mm)	Length (cm)
CETTA INF/N1.623CRV	PRIMO INF/N1.623CRV	1.6	23

Infusion Handle Non Valve

Comfortable handle to link tubing to needle, lack of valve reduces leakage. Single shaft allows infiltration and suction.

Product Code

CETTA INF/HAN W/O

Product Code

PRIMO INF/HAN W/O



Infusion Handle with Valve

This infusion handle, with on/off button, offers maximum control for infusion during the tumescent technique.

Product Code

PRIMO INF/HAN





Single Spike Tubing

Ideal for use with the HK pump, but will fit all devices 15.5 feet in length.

Latex and DEHP free.

Product Code

THK-ITS-20



Double Spike Tubing

Effort and time saving, two bags of solution at one time. Really good for treating multiple areas that require different solutions without bag changing.

15.5 feet in length.

Latex and DEHP free.

Product Code

THK-ITD-20



Tube Organiser

As the name explains, a very simple and effective way to organise single/multiple tubing during procedures.

X 100 pieces.

Product Code

THK-TO-100



Foot Pedal

Sealed and waterproof, additional foot pedals can be used with both the HK infusion pumps or HK aspirator. Useful to have one pedal either side of the operating table.

Product Code

THK-KIP-AP-FP





I.V. Pole and Bracket

Neatly and safely makes your infiltration system easily mobile around the operating theatre. Extremely stable with 6 legs and a very low centre of gravity.

Product Code

THK-IVPB



Infiltration Syringe – 30ml, 50ml

Spring loaded with a two way Luer lock valve attached, self-fills on retraction making hand held infusion faster and easier.

May we suggest for those with smaller hands the 30cc for comfort

Product Codes

Syringe - 30ml - INF

Syringe - 50ml - INF



I.V. Water Bath

A constant warmth water bath offer the advantages of obtaining consistent and controlled I.V. solution temperatures.

Product Code

THK-CTWB



Pressure Cuff Infuser

Simply wrap around an I.V. bag and inflate the cuff to keep outward pressure on the fluid. 3 way stop cock to prevent air leaks.

Product Code

TIN- INF-CUFF



Powered Suction



HK Aspirator

Super powerful performance, minimal noise.

The HK aspirator is a built for purpose liposuction machine of exceptional quality.

These much loved American style aspirators disappeared from the European Market for a few years as the USA manufacturers felt no inclination to acquire CE marking for our "little European market".

We kept receiving enquiries and the good news is..... A new CE marked beast of a Liposuction machine!

Weak suction, or suction slow to build up, not only slows down the procedure time but can result in increased tissue trauma.

The HK aspirator offers the power and peak suction speed so vital to this procedure.

Please get in contact for full spec details.

Product Code

THK-ADP-111

CE



Aspirator Disposable Tubing

Vital to successful power suction. Thick walled non collapsible non slip tubing made from a remarkably strong but light Polymer mixture.

3 Metre lengths, Latex free.

Product Code

FU 813300



Cetta - Disposable one piece cannulae

High grade stainless steel shaft, for aspiration with powered suction.



Cetta Handle: POM

Mercedes

Three holes spaced evenly around the tip.

The most popular tip style world-wide.



30 Hole Harvest Cannula

Harvest fat for transfer in smaller parcels for easier re-injection, reduces blocking.



Lantern

Chinese lantern shaped tip, useful for fibrous areas and popular in Gynecomastia procedures.



The Scream

Two holes above a single hole all downward facing. A cannula lookalike of the famous Munch painting.



One hole, flat, spoon backed for areas closer to the skin; such as chins, knees and ankles.




One Hole

A single downward facing aperture.



King Frank

Created by Dr Frank Niessen as a Multifunctional and multi useable cannula, holes on one face and solid back.



Product Code	Diameter (mm)	Length (cm)
CETTA 325/30holeM	3	25

Product Code	Diameter (mm)	Length (cm)
CETTA 325LM	3	25
CETTA 430LM	4	30
CETTA 530LM	5	30

Product Code	Diameter (mm)	Length (cm)
CETTA 210TSM	2	10
CETTA 320TSM	3	20
CETTA 430TSM	4	30

Product Code	Diameter (mm)	Length (cm)
CETTA 210SM	2	10
CETTA 315SM	3	15

Product Code	Diameter (mm)	Length (cm)
CETTA 2150HM	2	15
CETTA 3200HM	3	20
CETTA 4300HM	4	30

Product Code	Diameter (mm)	Length (cm)
CETTA 330KF/M	3	30
CETTA 430KF/M	4	30



Primo - Re-usable one piece cannulae

High grade stainless steel shaft, for aspiration with powered suction.

Primo Handle: Aluminium



PRIMO RANGE		
Product Code	Diameter (mm)	Length (cm)
PRIMO 210MM	2	10
PRIMO 315MM	3	15
PRIMO 325MM	3	25
PRIMO 430MM	4	30
PRIMO 420MM	4	20
PRIMO 530MM	5	30

Product Code	Diameter (mm)	Length (cm)
PRIMO 325/30holeM	3	25

Product Code	Diameter (mm)	Length (cm)
PRIMO 325LM	3	25
PRIMO 430LM	4	30
PRIMO 530LM	5	30

Product Code	Diameter (mm)	Length (cm)
PRIMO 210TSM	2	10
PRIMO 320TSM	3	20
PRIMO 430TSM	4	30

Product Code	Diameter (mm)	Length (cm)
PRIMO 210SM	2	10
PRIMO 315SM	3	15

Product Code	Diameter (mm)	Length (cm)
PRIMO 2150HM	2	15
PRIMO 3200HM	3	20
PRIMO 4300HM	4	30
Product Code	Diameter (mm)	Length (cm)
PRIMO 330KF/M	3	30
PRIMO 430KF/M	4	30

Mercedes

Three holes spaced evenly around the tip. The most popular tip style world-wide.



30 Hole Harvest Cannula

Harvest fat for transfer in smaller parcels for easier re-injection, reduces blocking.



Lantern

Chinese lantern shaped tip, useful for fibrous areas and popular in Gynecomastia procedures.



The Scream

Two holes above a single hole all downward facing. A cannula lookalike of the famous Munch painting.

Submental

One hole, flat, spoon backed for areas closer to the skin; such as chins, knees and ankles.



A single downward facing aperture.



King Frank

Created by Dr Frank Niessen as a Multifunctional and multi useable cannula, holes on one face and solid back.







Disposable Filter

The HEPA "High Alpha" 0.3 micron aspirator filter is a true high flow filter. The media is 144 square inches of pleated H.E.P.A fibreglass material. Under use with high powered aspirators the medial will not tear or obstruct. The housing is moulded from an unbreakable virgin polycarbonate. At a DOP rating of 99.975%, the HFIFD will filter up to 225 litres of air flow per minute.

Product Code

Filters



Disposable Liner

The liner is composed of a drainage receptacle and cap in a one piece unit, thus eliminating leakage and allowing disposal with no open contact.

Product Code

Liners



Abbott Canister

The Abbott canister system is the most popular collection system used to date and is simple and easy to use. Reinforced container shell to stop implosion.

Product Code

Canister



Tubing to Luer Lock Adapter

Small disposable adapter that allows Luer Lock cannula to be attached to machine suction.

Product Code

T530602



Toomey to Machine Suction Adaptor

Turn your Toomey cannula into machine suction devices using this simple syringe adaptor.

Product Code

TS-MSA





Medela Dominant Flex

The Dominant Flex is Medela's latest innovation in high-vacuum suction. One-touch adaptable flow rates allow for a wider variety of procedures to be carried out with just one pump. Designed for use in hospitals, clinics and medical practices, the Dominant Flex is the optimal medical suction pump for all suctioning needs.

Product Code

TMED-DOMFLEX



Medela Overflow Filter

To protect the pump from overflow, to protect the environment from bacteria or to neutralise odour.

Product Code

077.0571



Medela Canister 2.5L

Reusable PC jars offer a clean, simple handling and efficient fluid management.

Product Code

077.0085



Medela Disposable Liner 2.5L

Disposable suction 2.5 litre liners and reusable PC jars offer a clean, simple handling and efficient fluid management.

Product Code

077.0086



Medela Disposable Liner 2.5L with Solidifier

Disposable suction 2.5 litre liners with solidifier and reusable PC jars offer a clean, simple handling and efficient fluid management

Product Code

077.0087



Toomey Harvest

Toomey Cannula

High grade Stainless steel shaft, for aspiration with 60cc Toomey syringes.





Primo Hub: Aluminium, Stainless Steel Shaft

Cetta Hub: Aluminium, Stainless Steel Shaft

Mercedes Tip

Three holes spaced evenly around the tip.

The most popular tip style world-wide.



20/30 Hole Harvest Cannula

Harvest fat for transfer in smaller parcels for easier re-injection, reduces blocking.



Submental

One hole, flat, spoon backed for areas closer to the skin; such as chins, knees and ankles.



Open Harvest

Guarded fully open tip, designed specifically for tissue harvest.



King Frank

Created by Dr Frank Niessen as a Multifunctional and multi useable cannula, holes on one face and solid back



CETTA RANGE	PRIMO RANGE		
Product Code	Product Code	Diameter (mm)	Length (cm)
CETTA 210MTS	PRIMO 210MTS	2	10
CETTA 315MTS	PRIMO 315MTS	3	15
CETTA 318MTS	PRIMO 318MTS	3	18
CETTA 325MTS	PRIMO 325MTS	3	25
CETTA 430MTS	PRIMO 430MTS	4	30

Product Code	Product Code	Diameter (mm)	Length (cm)
CETTA 2.515/20holeTS	PRIMO2.515/20holeTS	2.5	15
CETTA 320/30holeTS	PRIMO320/30holeTS	3	20

Product Code	Product Code	Diameter (mm)	Length (cm)
CETTA 210SMTS	PRIMO 210SMTS	2	10
CETTA 315SMTS	PRIMO 315SMTS	3	15

Product Code	Product Code	Diameter (mm)	Length (cm)	
CETTA 318HARV/TS	PRIMO 318HARV/TS	3	18	

Product Code	Product Code	Diameter (mm)	Length (cm)
CETTA 330KF/TS	PRIMO 330KF/TS	3	30
CETTA 430KF/TS	PRIMO 430KF/TS	4	30





Toomey Syringe

60cc Toomey tip syringe – blister packed.

Product Code

TOOMEY-60CC

Toomey Syringe Locks



Snapper

60cc syringe lock - single arm. The simplest pull back one stage lock.

Product Code	Product Code
CETTA SNAP60	PRIMO SNAP60



Double Arm

60cc syringe lock - double arm. Base fits around the base of the syringe it helps to avoid droppage.

Product Code	Product Code
CETTA DOUBSNAP60	PRIMO DOUBSNAP60



2 Step Snap

60cc syringe lock - 2 step. Build pressure & suction in two stages.

Product Code	Product Code
CETTA2STEP	PRIMO2STEP



Square Back

60cc square syringe lock. Locks at any point along the syringe with a twist.

Product Code

PRIMOSQUARE60





Toomey to Luer Lock Adaptor

Transfer aspirate from 60cc Toomey to any size luer lock syringe.

Product Code	Product Code
CETTA MIC/LL/S	PRIMO MIC/LL/S



Toomey to Luer Lock Cannula Adaptor

Toomey syringe to luer lock needle adaptor

Product Code	Product Code	
CETTA MIC/LL/N	PRIMO MIC/LL/N	



Toomey to Toomey Fat Filter Hub & Syringe Stand

Stand includes 5 sets of 'O' rings and filters

Product Code	
PRIMO-GRAMS-FH	Filter hub & syringe stand
CETTA-GRAMSFHF5	Pack of 5 spare filters



Toomey Syringe Caps

Syringe caps, pack of 30.

Product Code	
PURE-SYR-CAP	



Toomey Syringe Rack 60cc

Autoclavable steel, holds 12 x 60cc syringes.

Product Code

PRIMO MIC/60CCR



Luer Lock Harvest

Luer Lock Cannula

High grade Stainless steel shaft, for aspiration with luer lock syringes.





Primo Luer Lock Hub, Stainless Steel, Shaft Stainless Steel





CETTA RANGE	PRIMO RANGE		
Product Code	Product Code	Diameter (mm)	Length (cm)
CETTA 215M/LL	PRIMO 215M/LL	2	15
CETTA 315M/LL	PRIMO 315M/LL	3	15
CETTA 318M/LL	PRIMO 318M/LL	3	18
CETTA 322M/LL	PRIMO 322M/LL	3	22
CETTA 325M/LL	PRIMO 325M/LL	3	25

Product Code	oduct Code Product Code		Length (cm)
CETTA 315/HARV/LL	PRIMO 315/HARV/LL	3	15
CETTA 322/HARV/LL	PRIMO 322/HARV/LL	3	22

Product Code	Product Code	Diameter (mm)	Length (cm)
CETTA 315S/LL	PRIMO 315S/LL	3	15
CETTA 322S/LL	PRIMO 322S/LL	3	22

Product Code	Product Code	Diameter (mm)	Length (cm)
CETTA 320KF/LL	PRIMO 320KF/LL	3	20
CETTA 330KF/LL	PRIMO 330KF/LL	3	30

1)
)

Mercedes Tip

Three holes spaced evenly around the tip.

The most popular tip style world-wide.



Open Harvest

Guarded fully open tip, designed specifically for tissue harvest.



Submental

One hole, flat, spoon backed for areas closer to the skin; such as chins, knees and ankles.



Created by Dr Frank Niessen as a Multifunctional and multi useable cannula, holes on one face and solid back.



20/30 Hole Harvest Cannula

Harvest fat for transfer in smaller parcels for easier re-injection, reduces blocking.





Cell Select

The ingenious Cell Select micro harvest cannula range is specifically designed to capture smaller aliquots of healthy fat cells via numerous 1 mm micro-ports. The harvest of smaller cell packages

of healthy fat cells via numerous 1mm micro-ports. The harvest of smaller cell packages results in easier, smoother re-injection through finer injection cannula. The spacing and increased surface area of the apertures per cannula tip results in an effective and constant collection speed, and greatly reduces blocking.

Product Code	Product Code	Holes	Diameter(mm)	Length(cm)	
CETTA 2515G	PRIMO 2515G	21	2.5	15	
CETTA 320G	PRIMO 320G	30	3	20	
LEVEL 2 - N	/ledium Rasp				
Product Code	Product Code	Holes	Diameter(mm)	Length(cm)	9 0 0
CETTA 2515M	PRIMO 2515M	21	2.5	15	
CETTA 320M	PRIMO 320M	30	3	20	
LEVEL 3 - A	Aggressive Ra	asp			
Product Code	Product Code	Holes	Diameter(mm)	Length(cm)	
CETTA 2515A	PRIMO 2515A	21	2.5	15	
CETTA 320A	PRIMO 320A	30	3	20	

Syringe Locks

Syringe lock - single arm

Product Code	Product Code	Syringe Size(cc)
CETTA SNAP10	PRIMO SNAP10	10cc
CETTA SNAP20	PRIMO SNAP20	20cc



Square syringe locks

Product Code	Syringe Size(cc)
PRIMOSOUARE10	10cc







VacLok Self Locking Syringe

Disposable self-locking syringe, locks at 5cc increments along the barrel.

Product Code

VAC-LOK 30cc



Luer Lock Syringes

Beckton & Dickinson syringe, luer lock, central nozel & latex free. Hand held harvest and injection.

Product Code

Syringe-10ml-BD

Syringe-20ml-BD



Luer Lock Syringe Cap

Soft, flexible cap to close syringe securely.

Product Code

Pure-Syr-Cap



Syringe Rack 10cc

Autoclavable steel rack, holds up to 60 x 10cc syringes

Product Code

MIC/10ccR



Fat Transfer

Re-Injection Cannula

High grade Stainless Steel for re-injection using Luer lock syringes



CETTA RANGE	PRIMO RANGE			
Product Code	Product Code	Diameter (mm)	Length (cm)	
CETTA MIC16/A7	PRIMO MIC16/A7	16G	7	
CETTA MIC16/A9	PRIMO MIC16/A9	16G	9	
CETTA MIC16/B7	PRIMO MIC16/B7	16G	7	
CETTA MIC16/B9	PRIMO MIC16/B9	16G	9	
CETTA MIC16/B12	PRIMO MIC16/B12	16G	12	
CETTA MIC16/B18	PRIMO MIC16/B18	16G	18	
CETTA MIC16/C7	PRIMO MIC16/C7	16G	7	
CETTA MIC16/C9	PRIMO MIC16/C9	16G	9	
CETTA MIC18/A9	PRIMO MIC18/A9	18G	9	
CETTA MIC18/B9	PRIMO MIC18/B9	18G	9	



Curved Cannula

Product Code	Product Code	Diameter (mm)	Length (cm)
CETTA CURV/18B7	PRIMO CURV/18B7	18G	7
CETTA CURV/16B12	PRIMO CURV/16B12	16G	12
CETTA CURV/16B15	PRIMO CURV/16B15	16G	15

Ultra Fine

Ultra fine facial needle with unblocking Stylet







Curved obtuse tip for accurate placement of micro droplets

Product Code	Product Code	Diameter (mm)	Length (cm)
CETTA 1285MD	PRIMO 1285MD	1.2	8.5





Dissectors

V-Dissector

Undermining fibrous tissue, widens at the tip and ideal for breast surgery.

CETTA RANGE	PRIMO RANGE		
Product Code	Product Code	Diameter (mm)	Length (cm)
CETTA VDISS167	PRIMO VDISS167	16G	7
CETTA VDISS169	PRIMO VDISS169	16G	9



W-Dissector

For undermining small areas of fibrous tissue, less aggressive and finer than V Dissector.

Product Code	Product Code	Diameter (mm)	Length (cm)
CETTA WDISS	PRIMO WDISS	16G	8.5



Hay Tunnel Dissector

For blunt dissection tunnelling when a V Dissector may be considered aggressive i.e. around implants. Slight upward curve to tip in order to inhibit direct forceful contact, comfortable handle with orientation marking.

Product Code	Product Code	Diameter (mm)	Length (cm)
CETTA HTD	PRIMO HTD	2	12
32.0			The second second second



Injection Guns

Cell Brush 3cc/10cc

Precise micro droplets delivery, removes pressure allowing even, controlled delivery. 10ml syringe delivers 0.50ml with each brush of the operators thumb. 3ml = 0.15 per brush.

Product Code	Capacity (cc)
CELLBRUSH3	3
CELLBRUSH10	10

The Golden Gun 10cc

Thumb controller injection from 10cc BD syringe, eases the physical pressures of re-injection and controls accidental over filling.



Trigger Injection Gun 10cc/5cc

Trigger mechanism controlled injection from 10cc BD syringe or 5cc BD syringe (adaptor included) eases the physical pressures of re-injection and controls accidental over filling.









Micro/Nano Filters

Want to reduce fat molecule size or filter fat for re-injection?



Luer Lock to Luer Lock Transfers - Standard Hole

Simply use as a standard transfer adaptor tool between luer lock syringes.

Or for finer anaerobic re-injectable tissue start the process by working the graft through the adaptor from syringe to syringe around 20 times.

Product Code	Product Code
CETTA MIC/LL/LL	PRIMO MIC/LL/LL

Finer graft required, then:



Luer Lock to Luer Lock Micro Adaptors

Decrease fat size further in stages using these specially drilled multi-holed adaptors for a graft consistency more easily injectable via fine facial needles.

Primo only

Product Code

PRIMO MFC/LL

Then finer still:



Nano Fat Mesh Filter

This is designed to be used for very fine injectable grafts.

Re-usable outer case with sterile CETTA disposable mesh sets.

Product Code

PRIMO NFC/LL Adaptor

CETTA NFDF Mesh x 5 per pack



Autologous Fat Filters

Popular and easy to use Container Systems



Aquavage is available in 1 or 2 litre containers, this is a sterile, single-use, closed system for high volume fat harvesting sitting as a link between patient and suction pump.

The Aquavage container uses fluid to float and separate adipose tissue, this allows the removal first of the fluid portion and then the fat cells from the container via a sterile and closed method. Fat is decanted from the bottom of the Aquavage system and can be injected straight back to the patient – no need for a centrifuge reducing any contaminants.

Product Codes

AQUAVAGE - 1200

AQUAVAGE - 2000

AQUAVAGE- STAND

Filtron

We supply 100cc, 250cc 500cc.

Standard tumescent solution is infused. The Filtron unit is placed in line between the harvest cannula and the waste canister/suction source. Large volume units may be attached by the universal bracket to a sterile-draped canister stand. Small volume units may be placed in the sterile field. Suction pressure is reduced to 15"Hg and fat is drawn into the mesh filter within the canister. Waste and non-viable cells are pulled through the filter and are suctioned out to the waste canister via the removal tube. When the desired volume is collected, the FILTRON unit is then inverted and the transfer tube attached. The quality fat harvest is then drawn out directly into a luer lock syringe and is ready to re-inject. Sterile, fast, effective, quality harvest without exposure to ambient air.



Product Codes

FILTRON/ 100

FILTRON/250

FILTRON/500



Product Code

CENTRIFUGE

The MPW 54 Centrifuge 50/60Hz- Separation of fat or bloods for PRP

Speed: 3500, 5800 rpm or 1000^* , 3500^* rpm RCF: 1137, $3122 \times g$ or 90^* , $1137^* \times g$

Dimensions: 180 x 220 x 270mm Weight of centrifuge 230 V: 4kg

Features

- LCD display
- 2 centrifuging speeds
- 1 acceleration/deceleration curve
- Easy operation
- Centrifuging time adjustment
- Continuous operation mode
- Small size
- Autoclavable barrel and tubes.

Safety

- Lid locking during rotor run
- Emergency lid lock release



PURE 3 RAFT

A closed sterile system for the selective filtration of autologous grafts

The Puregraft filtration bags offer standardisation of result with ease of use.

A very simple and ingenious system that offers superior separation of graft and allows the gentle preservation of the regenerative properties of fat.

- 1) Simply add syringe harvested lipo aspirate
- 2) Wash tissue inside the bag with added lactated ringers
- 3) Drain away filtered waste, including free lipids and contaminants
- 4) Collect injection ready graft from the bag with your choice of luer lock syringe



RAISING
THE STANDARD
FOR FAT
GRAFTING

Product Codes

Pure-50/PURE Grafts 0 - 100cc

Pure-250/PURE Grafts 50cc - 400cc

Pure-850/PURE Grafts 200cc - 1200cc




Cost effective, choose the size of bag that fits the requirements

Puregraft provided in partnership with









Fully closed tubing and adaptor system for collection and injection of graft without syringes.

Red - One way harvest

Yellow - One way injection

Product Code

PURE-FLS

Tubing and adaptors

Easel & Slider Set

For 250cc, 850cc

Allows easy positioning and superior gravity separation during procedure

Product Code

PURE-550/IS



Let us make life easy for you, we will create a unique kit for each surgeon/procedure and code it accordingly for problem free ordering.

i.e. Mr Smith / Small Volume Kit

- PG50 50cc Puregraft Filter
- VacLok self-locking syringe
- FatLok Adaptor / tubing system
- Cell Select Infusion / harvest cannula
- 15FN Facial injection cannula
- 5 x 1ml BD luer lock syringes





Compression Garments

Arm garment

CLG400 - Medium Sleeve Style Arm Cosmetic Surgery Compression

Applicable Procedures: Arm Lift, Arm Liposuction

Features: Designed for surgical procedures of the arm, with Velcro closure in the back and front. The sleeve length is slightly below the elbow.

- Provides optimal amount of compression to help improve the blood circulation, minimizes swelling after the procedure, and speeds up the healing process.
- Provides support to surgical areas for more comfort and helps the skin fit better to its new contours.
- Flushes the body out of potentially harmful fluids, allows patient to return to daily routines sooner.
- The garment can be worn with a full body garment for added upper body compression and has adjustable Velcro closures for more comfort and ease of application.
- Appropriate for liposuction and arm procedures (lift, reduction, brachioplasty), as well as for arm procedures following gastric bypass or bariatric surgery.
- Designed with Velcro closure on the back and front for extra support and compression around the shoulders and back area.
- Constructed with a high performance Comprexx TecSheen fabric that is long lasting and retains the garments shape and compression even under extended wear and washing.





Size Guide

	S	М	L	XL
ARM	10.25-11 IN	11-12 IN	12-13 IN	13-14.25 IN
	26-28 CM	28-30.5 CM	30.5-33 CM	33-36 CM

Facial garment

Size Guide

	S	М	L	XL
NECK (Unisex)	13-15 IN	15-17 IN	17-19 IN	19-21 IN
	33-38 CM	38-43 CM	43-48 CM	48-53 CM

Applicable Procedures: Browlift, Face Lift, Neck Lift, Liposuction (face, chin, neck), Chin & Jaw Augmentation, Genioplasty, Facial Implants, Anti-Snoring, Controlling Dry Mouth, Trauma Compression, Maxillofacial Surgery, Platysmaplasty, Oral Surgery, Coronal Lift, Submentoplasty, Submandibular Lipectomy, Mastoid Surgery, Cervicofacial Liposuction.

Features: UNISEX smooth compression facial garment designed chin strap with Velcro closure on top and medium neck support. This compression facial garment acts like a "dri-fit" breathable, latex free and antimicrobial which prevents odour and bacteria build up.

- Provides compression in the face, chin, neck areas and is appropriate for a face lift, neck lift, and liposuction.
- Improves blood circulation, minimises swelling after the procedure, flushes the body out of potentially harmful fluids, accelerates the healing process, and allows the patient to return to daily routines sooner.
- Ideal for use following a Face Lift, Neck Lift, Liposuction (face, chin, neck), Chin & Jaw Augmentation, Genioplasty, or Facial Implants surgery.
- Also provides support to surgical areas for more comfort and helps the skin fit better to it new contours.
- Features hook engage-able material on the outside enables you to adjust the right amount of pressure or support you feel most comfortable with.
- Easy pull on adjustable Velcro neck closures with ear open.
- Designed with a non-slip dome and a perfectly placed chin cup for continuous compression and support.
- Made of woven elastic anti-roll material for more durable compression and performance and has a soft plush side that feels comfortable against the skin. Moreover, the garment is made of pattern designs that target maximum support and compression.





CLG100 - Chin Strap Support Compression Garment W/ Medium Neck Support




Abdominal binders (unisex)

CLG1350 - 12" Abdominal Plastic Surgery 4 Split Panel Binder

Applicable Procedures: Liposuction of the abdomen and flanks, Tummy Tuck, Lipoplasty, Hernia Surgery, Back

Features: Designed with a Two panel Velcro closure.

- Made up of two 6" SPLIT panels (12" wide) for better flexibility which helps the skin fit better to it's new contours.
- Typically used by people with longer torsos following liposuction of the abdomen and flanks, tummy tuck or hernia surgery.



- Also provides support to surgical areas for more comfort and helps the skin fit better to it new contours.
- Constructed with woven elastic anti-roll material for more durable compression and performance.
- Features a soft plush side that feels comfortable against the
- Garments have a hook engage-able material on the outside enables you to adjust the right amount of pressure or support you feel most comfortable with. It also allows the Velcro to attach anywhere on the binder for precise fit and even compression.
- Gives a soft & flexible Velcro hook for a tight snug fit that lasts.







	R	XL
WAIST	23-45 IN 58.5-114.5 CM	45-62 IN 114.5-157.5 CM



CLG1250 - 12" Abdominal Plastic Surgery 4 Panel One Piece Binder

Applicable Procedures: Liposuction of the abdomen and flanks, Tummy Tuck, Lipoplasty, Hernia Surgery, Back Support.

Features: Designed with a Four panel Velcro closure.

- Made up of four 3" panels (12" wide) for better flexibility which helps the skin fit better to it's new contours.
- Typically used by people with longer torsos following liposuction of the abdomen and flanks, tummy tuck or hernia surgery.



CLG1300 - 9" Abdominal Plastic Surgery 3 Split Panel Binder

Applicable Procedures: Liposuction of the abdomen and flanks, Tummy Tuck, Lipoplasty, Hernia Surgery, Back Support.

Features: Designed with a three panel Velcro closure, three

- Made up of three 3" SPLIT panels (9") for better flexibility which helps the skin fit better to it's new
- Typically used by women or people with shorter torsos following liposuction of the abdomen and flanks, tummy tuck or hernia surgery.







CLG1200 - 9" Abdominal Plastic Surgery One Piece Binder

Applicable Procedures: Liposuction of the abdomen and flanks, Tummy Tuck, Lipoplasty, Hernia Surgery, Back Support.

Features: Designed with a three panel Velcro closure.

- Made up of three 3" panels (9") for better flexibility which helps the skin fit better to it's new contours.
- Typically used by women or people with shorter torsos following liposuction of the abdomen and flanks, tummy tuck or hernia surgery.





Male abdominal

CLG210 - Male High-Waist Abdominal Cosmetic Surgery Compression Brief with Zippers

Applicable Procedures: Abdominal Liposuction

CLG215 - Male Above Knee High-Waist Abdominal Cosmetic Surgery Compression Girdle with Zippers

Applicable Procedures: Abdominal & Thigh Liposuction

Features: Designed with a 3" high waist elastic band front zipper and pouch opening in crotch

- Provides optimal amount of compression to help improve the blood circulation, minimises swelling after the procedure, and speeds up the healing process.
- Provides support to surgical areas for more comfort and helps the skin fit better to its new contours.
- Flushes the body out of potentially harmful fluids, allows patient to return to daily routines sooner.
- The anti-microbial treatment also helps retard growth of bacteria or fungus that can cause odor skin irritation, mildew or fabric degradation.
- Features a plush 3" elastic band and a pouch opening on the crotch area.
- Provides comfort, support and compression and is designed with a front zipper.
- Constructed with a high performance Comprexx Tecsheen fabric that is long lasting and retains the garments shape and compression even under extended wear and washing.



Size Guide

Olze Guide					
	S	М	L	XL	XXL
CHEST	34-38 IN	38-42 IN	42-46 IN	46-50 IN	50-54 IN
	86.5-96.5 CM	96.5-106.5 CM	106.5-117 CM	117-127 CM	127-137 CM
WAIST	28-32 IN	32-36 IN	36-40 IN	40-44 IN	44-48 IN
	71-81 CM	81-91.5 CM	91.5-101.5 CM	101.5-112 CM	112-122 CM



CLG310 - Male Abdominal Cosmetic Surgery Compression Vest CLG315 - Male Abdominal Cosmetic Surgery Compression Vest With Sleeves

Applicable Procedures: Abdominal and Chest Liposuction, Augmentation and Reconstruction

Features: Designed with a zipper in front, Velcro shoulder straps, and without sleeves

- This compression vest is designed ideal for use following Plastic Surgery, liposuction of the abdomen and flanks and is typically worn during Stage 1 of recovery (initial 1-2 weeks of post-op).
- Provides support to surgical areas for more comfort and helps the skin fit better to its new contours.
- Provides optimal amount of compression to help improve the blood circulation, minimises swelling after the procedure, and speeds up the healing process.
- Flushes the body out of potentially harmful fluids, allows patient to return to daily routines sooner.
- The anti-microbial treatment also helps retard growth of bacteria or fungus that can cause odor skin irritation, mildew or fabric degradation.
- Features a soft & flexible Velcro shoulder straps.
- Provides comfort, support and compression and is designed with a front zipper.
- Constructed with a high performance Comprexx Tecsheen fabric that is long lasting and retains the garments shape and compression even under extended wear and washing.



Breast support

CLG200 - Breast Surgery Support Bra W/1" Elastic Band

Applicable Procedures: Breast Augmentation (Breast Implants, Transaxillary Submuscular Breast Augmentation, Trans Umbilical Breast Augmentation), Breast Lift (Mastopexy), Breast Reduction, Breast Reconstruction, Latisimus Dorsi procedure, Mastectomy, TRAM Flap Procedure.









- Ideal for use following plastic surgery and any breast procedure, this surgical bra provides optimal amount of compression to help improve the blood circulation, minimises swelling after the procedure, and speeds up the healing process.
- Provides support to surgical areas for more comfort and helps the skin fit better to its new contours.
- Flushes the body out of potentially harmful fluids, allows patient to return to daily routines sooner.
- The anti-microbial treatment also helps retard growth of bacteria or fungus that can cause odor skin irritation, mildew or fabric
- It is constructed with a high performance comprexx tecsheen fabric.
- The Comprexx Tecsheen fabrics are long lasting and retain their shape and compression even under extended wear and washing.
- Features a plush 1" elastic band with adjustable hook & eye shoulder straps.
- Also features a front hook & eye opening.

Size Guide

	S	М	L	XL	XXL
UNDER BUST	24-28 IN	28-32 IN	32-36 IN	36-40 IN	40-44 IN
	61-71 CM	71-81 CM	81-91 CM	91-101.5 CM	101.5 - 112 CM

CLG208 - Seamless Bra

Applicable Procedures: Fat Transfer, Breast Augmentation (Breast Implants, Transaxillary Submuscular Breast Augmentation, Trans Umbilical Breast Augmentation), Breast Lift (Mastopexy), Breast Reduction, Breast Reconstruction, Latisimus Dorsi procedure, Mastectomy, TRAM Flap Procedure.

Features: Designed with no cup seam for maximum comfort.









CLG207 - 3" Wide Breast Surgery Implant Stabilizer Band

Applicable Procedures: Breast Augmentation (Breast Implants, Transaxillary Submuscular Breast Augmentation, Trans Umbilical Breast Augmentation), Stabilising & Positioning Implants.

Features: Designed with Velcro closure, one size fits all.

One Size Fits All



- Ideal for use following a breast augmentation.
- Made of mono filament segmented breathable elastic.
- Can be used alone or with a bra to enhance post-operative results of breast augmentation procedures.
- Features adjustable Velcro closures that enable customised band positioning and compression to either breast.
- Helps stabilise the breast implant, keep it in its correct position, and prevent the implant from rising.
- The Adjustable Breast Surgery Implant Stabiliser Band is universal in size.



Girdles & body suits

Size Guide

	S	М	L	XL	XXL
WAIST	25-27 IN	27-29.5 IN	29.5-32.5 IN	32.5-36 IN	36-40 IN
	64-69 CM	69-75 CM	75-83 CM	83-97 CM	97-102 CM
HIP	35.5-37.5 IN	37.5-40 IN	40-43 IN	43-47 IN	47-52 IN
	90-95 CM	95-102 CM	102-109 CM	109-119 CM	119-132 CM
THIGH	19.5-21 IN	21-22.5 IN	22.5-24 IN	24-25.5 IN	25.5-27 IN
	50-53 CM	53-57 CM	57-61 CM	61-65 CM	65-69 CM



High Waist Abdominal Compression Girdle

Applicable Procedures: Liposuction in the Abdomen, Buttock Lifts, Tummy Tuck

Features: Designed with a zipper on the left side, hook and eye crotch flap with a non roll 3" waist elastic band



CLG500 - Above the

Knee











Ideal for use following Plastic Surgery, liposuction of the, thighs, hips, buttocks, and

- Worn during Stage 1 of recovery (initial 1-2 weeks of post-op).
- Provides optimal amount of compression to help improve the blood circulation, minimises swelling after the procedure, and speeds up the healing process.
- Provides support to surgical areas for more comfort and helps the skin fit better to its new contours.
- Flushes the body out of potentially harmful fluids, allows patient to return to daily routines
- The anti-microbial treatment also helps retard growth of bacteria or fungus that can cause odor skin irritation, mildew or fabric degradation.
- Features a 3" plush waist elastic band, a front double panel for extra support and compression and a padded hook & eye crotch flap.
- The Comprexx Tecsheen fabrics are long lasting and retain their shape and compression even under extended wear and washing.



CLG900 - Ankle Lenath

CLG300 - Panty Length

Compression Body Suit With Bra

Applicable Procedures: Lower Body Lift, Buttocks Augmentation, Implants, Lift (Gluteplasty), Liposuction (Upper and Lower Abdomen, Back, Flank, Hip, Thigh), Lipoplasty, Liposculpture, Lower Body Procedures, Plastic surgery procedures following Gastric Bypass Surgery, Bariatric Surgery, and Stomach Stapling, Tummy Tuck (Abdominoplasty).

Features: Designed with a built-in surgical compression bra, zipper on the left side, high back compression support of the back and a hook and eye crotch flap for ease. Smooth compression garment that is like a "dri-fit" breathable, latex free and antimicrobial which prevents odor and bacteria build up. Recommended use during the first 1-2 weeks post surgery.















Measuring guides

MEASURING INSTRUCTIONS:

Hold measuring tape against skin while measuring

Do not pull tape tightly against the skin

Do not let measuring tape snag

(If you are between sizes, you may prefer the larger size).

Sleeves

	S	М	L	XL
ARM	10.25-11 IN	11-12 IN	12-13 IN	13-14.25 IN
	26-28 CM	28-30.5 CM	30.5-33 CM	33-36 CM



Facial Garments

	S	М	L	XL
NECK (Unisex)	13-15 IN	15-17 IN	17-19 IN	19-21 IN
	33-38 CM	38-43 CM	43-48 CM	48-53 CM



Male Abdominal & Male Brief

	S	М	L	XL	XXL
CHEST	34-38 IN	38-42 IN	42-46 IN	46-50 IN	50-54 IN
	86.5-96.5 CM	96.5-106.5 CM	106.5-117 CM	117-127 CM	127-137 CM
WAIST	28-32 IN	32-36 IN	36-40 IN	40-44 IN	44-48 IN
	71-81 CM	81-91.5 CM	91.5-101.5 CM	101.5-112 CM	112-122 CM





Garment colours. Please note, the images shown may not represent the true colour of the available garment style, please enquire.

Abdominal Binders

	R	XL
WAIST	23-45 IN 58.5-114.5 CM	45-62 IN 114.5-157.5 CM



Breast Support

	S	М	L	XL	XXL
UNDER BUST	24-28 IN 61-71 CM	28-32 IN 71-81 CM	32-36 IN 81-91 CM	36-40 IN 91-101.5 CM	40-44 IN 101.5-112 CM



Girdles & Body Suits

	S	М	L	XL	XXL
WAIST	25-27 IN	27-29.5 IN	29.5-32.5 IN	32.5-36 IN	36-40 IN
	64-69 CM	69-75 CM	75-83 CM	83-97 CM	97-102 CM
HIP	35.5-37.5 IN	37.5-40 IN	40-43 IN	43-47 IN	47-52 IN
	90-95 CM	95-102 CM	102-109 CM	109-119 CM	119-132 CM
THIGH	19.5-21 IN	21-22.5 IN	22.5-24 IN	24-25.5 IN	25.5-27 IN
	50-53 CM	53-57 CM	57-61 CM	61-65 CM	65-69 CM





Garment Care











Ancillary Equipment



Punches

Punch incisions for healing without suturing. Also useful for acne scar treatment.

Product Code

BP15	1.5mm x 20 per pack
BP2	2mm x 20 per pack



PRP Kit

A very simple little disposable kit allowing the closed preparation of concentrated PRP for re-injection with or without fat graft.

- 1 Vacuum Tube 10ml
- 1 Sleeve Filter
- 1 Blunt Needle 18x100mm
- 1 10ml Syringe (to be used with the Blunt Needle)
- 1 VACU10® PSV Set
- 1 1ml luer lock Syringe
- 1 Hypodermic Needle 30G



Cannula Cleaning Brushes

Product Codes

Product Code

PRP ECO KIT

TS-BR2	Cannula cleaning brush, 2mm x 42cm
TS-BR3	Cannula cleaning brush, 3mm x 48cm
TS-BR4	Cannula cleaning brush, 4mm x 48cm
TS-BR5	Cannula cleaning brush, 5mm x 48cm
TS-BR6	Cannula cleaning brush, 6mm x 48cm
TS-CCS1630	Cannula cleaning stylet,16ga x 30cm
TS-CCS1830	Cannula cleaning stylet, 18ga x 30cm



Syringes

Product Codes

Syringe-1ml-BD
Syringe-3ml-BD
Syringe-5ml-BD
Syringe-10ml-BD

Product Codes

Syringe-20ml-BD
Syringe-50ml-BD
Syringe-Softjet 10ml - CELLBRUSH
Syringe-Softjet 3ml - CELLBRUSH





Midline Position Pillow

The Thigh MidlineTM is a wedge-shaped surgical positioning pillow that helps approximate the anatomic area with the patient in the lateral decubitus position. An essential tool for outer thigh liposuction procedures.

The Midline Sterile Plastic Covers are used in conjunction with the Midline positioning pillow and are single-use.

Product Code

MPPIL Thigh Midline Pillow x 1

MPPIL-SC Sterile Plastic Cover x 25 per pack



Lipo Positioning Pillow

The Thigh Aside™ pillow helps position and maintains the targeted thigh area in the ideal anatomic position during liposuction of inner thighs and inner knees.

The Aside Sterile Plastic Covers are single-use and sterile plastic covers used in conjunction with the Aside positioning pillow.

Product Code

TAPIL Thigh Aside Pillow x 1

TAPIL-SC Sterile Plastic cover x 25 per pack



Photo Panty

Disposable latex free fabric underwear.

Standardise and improve pre and post-operative photographs.

Provide comfort and dignity during examination for a very low cost.

Product Code

PIC-EZE



Nourisil Scar Gel

NourisilTM MD contains a unique blend of five silicones and Vitamin E to improve the appearance of scars and maintain the skin's moisture balance. NourisilTM MD has been designed for the treatment and prevention of hypertrophic and keloid scars which result from medical and cosmetic surgery, wounds, trauma or burns.

Product Code

ScarGel-30ml

Expertise • Experience • Guidance

Triangle Surgical is a specialist medical distribution company in the sector of Cosmetics, Aesthetics and Reconstruction. Our focus is on designing, sourcing and marketing superior products that offer real advantage and advancement to the doctor.

We are a long established family run company with significant combined experience in this field and we take great pleasure and pride in our team's outstanding reputation for quality of service and in-depth product knowledge.

www.trianglesurgical.uk

Master Indemnity Registration number: DHMIA/1773/17 UK

Scotland: MIA 1773/17 Wales: 1773/17 NHSW

EC DIRECTIVE 93/42/EEC+A: 2007/47/EC
EC DIRECTIVE 93/42/EEC ANNEX V ARTICLE 3

EN ISO 134485:2003 + AC2009

EN ISO 9001:2008 EN ISO 1137-1:2006

EN ISO 13485:2016+AC 2016

Pending: EN ISO 13485: 2016 MDR

6-10 Church Farm Courtyard High Street Chalfont St Giles Buckinghamshire HP8 4QH

Tel: +44 (0) 1494 876990 Fax: +44 (0) 1494 876991 Email: info@trianglesurgical.uk

