

Our most  
**versatile**  
LED device

Discover  
**COMPACT LITE**



 DERMALUX

# Precise, proven and versatile



Introducing the **Dermalux Compact Lite**, a **revolutionary** tabletop professional LED device, powered by our **unique technology** for **consistent and unrivalled results**.

Compact Lite delivers the **most clinically proven wavelengths**, backed by science to target a wide range of conditions **for face and body**.

The adjustable LED panel design enables **versatile positioning** for use with a treatment couch or chair.

## Life-improving innovations

**Improving people's lives** drives us to innovate. Our pioneering devices **harness our bodies' incredible powers to heal and repair** by targeting the deepest skin cells **to treat the true root cause of your patients' concerns**. The ability to transform people's skin health, enhance their wellbeing, quality of life even, is our motivation to keep **leading the way**.

## Compact Lite at a glance

- 1 The most clinically proven wavelengths: Blue 415nm, Red 633nm, Near Infrared 830nm
- 2 Next generation ceramic LED technology for maximum power and unrivalled results
- 3 Upgradeable software to future proof your investment
- 4 Advanced preset, post treatment and custom protocols
- 5 Use with a treatment bed or chair
- 6 Includes training, ongoing education and customer support

## Technical specification

Wavelengths (+/- 2NM)	<b>415NM, 633NM, 830NM</b>
Classification	<b>Non-Medical Device</b>
Treatment Time (Mins)	<b>1 - 20</b>
Dose Range (Joules)	<b>7 - 140</b>
Intensity (M/W)	<b>7 - 116</b>
Device Dimensions (MM)	<b>H: 700 W: 550 D: 360</b>
Device Weight (KG)	<b>32</b>





Want to learn more?  
Visit our website or  
speak to our team



dermalux.com

@dermaluxled



211 Europa Blvd, Westbrook, Warrington WA5 7TN



 DERMALUX